

IHE Integration Statement

| Vendor | Product Name | Version | Date |
|---------------------------|---|---------|-----------|
| Siemens Medical Solutions | SIENET Magic Integrated Radiology Suite | V42 | Nov. 2004 |

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

| Integration Profiles Implemented | Actors Implemented | Options Implemented |
|--|--|--|
| Scheduled Workflow (SWF) | Image Display Evidence Creator ^{*1)} Image Archive Image Manager PPS Manager ^{*2)} | no options defined none none none none |
| Patient Information Reconciliation (PIR) | Image Archive Image Manager PPS Manager ^{*2)} | none none none |
| Consistent Presentation of Images (CPI) | Print Composer ^{*1)} | No options defined |
| Access to Radiology Information (ARI) | Image Display | None |

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

| | | | |
|--------------|---|--|--|
| HL7 | n.a. | | |
| DICOM | http://www.siemens.com/dicom | | |

Links to general information on IHE

| | | |
|---|--|---|
| In North America: http://www.rsna.org/IHE | In Europe: http://www.ihe-europe.org | In Japan: http://www.jira-net.or.jp/ihe-j |
|---|--|---|

1) Component: Intelligent Postprocessing Workplace.

2) Includes Mitra Broker