

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	Somatom Huan Yue	VB10	09/2004

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query
		Broad Worklist Query
		Modality Group Case
	Image Display	No options defined in IHE
	Evidence Creator	none
Patient Information Reconciliation	Acquisition Modality	No options defined in IHE
Access to Radiology Information	Image Display	none
Consistent Presentation of Images	Print Composer	No options defined in IHE

**Internet address for vendor's IHE information:** <http://www.siemens.com/ihe>

### Links to Standards Conformance Statements for the Implementation

<b>HL7</b>	n.a.		
<b>DICOM</b>	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		

### Links to general information on IHE

In North America: <a href="http://www.rsna.org/IHE">http://www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>
---	--	---