

IHE Integration Statement

| Vendor | Product Name | Version | Date |
|---------------------------|----------------------|---------|---------|
| Siemens Medical Solutions | Somatom Sensation 10 | VB10 | 09/2004 |

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

| Integration Profiles Implemented | Actors Implemented | Options Implemented |
|------------------------------------|----------------------|------------------------------|
| Scheduled Workflow | Acquisition Modality | Patient Based Worklist Query |
| | | Broad Worklist Query |
| | | Modality Group Case |
| | Image Display | No options defined in IHE |
| | Evidence Creator | none |
| Patient Information Reconciliation | Acquisition Modality | No options defined in IHE |
| Access to Radiology Information | Image Display | none |
| Consistent Presentation of Images | Print Composer | No options defined in IHE |

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

| | | | |
|--------------|---|--|--|
| HL7 | n.a. | | |
| DICOM | http://www.siemens.com/dicom | | |

Links to general information on IHE

| | | |
|---|--|---|
| In North America: http://www.rsna.org/IHE | In Europe: http://www.ihe-europe.org | In Japan: http://www.jira-net.or.jp/ihe-j |
|---|--|---|