

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	Somatom Sensation 64	VB20	10/2005
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query	
		Broad Worklist Query	
		Modality Group Case	
	Image Display	No options defined in IHE	
	Evidence Creator	none	
Patient Information Reconciliation	Acquisition Modality	No options defined in IHE	
Access to Radiology Information	Image Display	none	
Consistent Presentation of Images	Print Composer	No options defined in IHE	
Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.rsna.org/IHE	In Europe: http://www.ihe-europe.org	In Japan: http://www.jira-net.or.jp/ihe-j	