

# IHE Integration Statement

Product Name: SC2000

Version: 2.0

Date: December 08, 2011

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	SC2000	2.0	Dec 08, 2011
<p>This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:</p>			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query	
		Patient Based Worklist Query	
Cardiology Evidence Documents	Evidence Creator	Echo Evidence Vascular Evidence	
Cardiology Evidence Documents	Acquisition Modality	Echo Evidence Vascular Evidence	
<p><b>Internet address for vendor's IHE information:</b> <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a></p>			
<p><b>Links to Standards Conformance Statements for the Implementation</b></p>			
<b>HL7</b>	n.a.		
<b>DICOM</b>	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
<p><b>Links to general information on IHE</b></p>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>		In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>