

IHE Integration Statement		Date	Jan, 2008
Vendor	Product Name	Version	
Siemens	Siemens X300	3.0	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Evidence Documents	Evidence Creator	OB-GYN Ultrasound Procedure Report (TID 5000) Adult Echocardiography Procedure Report (TID 5200) Vascular Ultrasound Procedure Report (TID 5100)	
Scheduled Workflow	Acquisition Modality	<ul style="list-style-type: none"> - Broad Worklist Query - Patient based Worklist Query 	
Echocardiography Workflow	Acquisition Modality	<ul style="list-style-type: none"> - Broad Worklist Query - Patient based Worklist Query 	
Patient information Reconciliation	Acquisition Modality	<ul style="list-style-type: none"> - None 	
Internet address for vendor's IHE information:		http://www.siemens.com/IHE	
Links to Standards Conformance Statements for the Implementation			
HL7	N/A		
DICOM	http://www.siemens.com/DICOM		
Links to general information on IHE			
In North America: www.ihe.Net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	