

IHE Integration Statement

Product Name:

syngo Dynamics Workplace

Version: 9.0

Date: December 2010

IHE Integration Statement		Date	01 Dec 2010
Vendor	Product Name		Version
Siemens Medical Solutions	syngo® Dynamics Workplace		9.0
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Displayable Reports (DRPT - CARD)	Report Creator	■	
Evidence Documents (ED - CARD)	Evidence Creator	<ul style="list-style-type: none"> ■ Echo Evidence <ul style="list-style-type: none"> ○ Vascular (5100) ○ Echocardiography (5200) 	
Retrieve ECG for Display (ECG - CARD)	Display	<ul style="list-style-type: none"> ■ Retrieve ECG List (CARD-5) ■ Retrieve ECG Document for Display (CARD-6) <ul style="list-style-type: none"> ○ Scalable Vector Graphics 	
Retrieve Information for Display (RID - ITI)	Display	■	
Internet address for vendor's IHE information:		http://www.siemens.com/IHE	
Links to Standards Conformance Statements for the Implementation			
HL7	N/A		
DICOM	http://www.siemens.com/DICOM		
Links to general information on IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	