

3T Unlimited




**MAGNETOM Trio**  
Easy and cost-efficient siting

# MAGNETOM Trio

## Easy and cost-efficient siting

True to its tradition of innovations, Siemens has continuously worked to facilitate the siting of MR systems. Through miniaturization of components at every level, the new generation of MAGNETOM® requires less space for the examination room. But what is even more important, less space is required for the electronic components as well.

MAGNETOM Trio stands for "3T Unlimited" and is the shortest 3T MR system with the largest (40 cm) FoV in all three axes – for whole-body applications.



### Features

- Min. area < 33 m<sup>2</sup> (< 357 sq. ft.)
- Min. room height clearance 2.73 m (8'11")
- Only 4 electronic cabinets, including 3 that do not require access from the sides or from the back
- Power connection value 95 kVA

### Benefits

- Reduced siting costs
- Low power requirements
- More space for other requirements

#### Please contact in USA:

Siemens Medical Solutions USA, Inc.  
51 Valley Stream Parkway  
Malvern, PA 19355  
Tel.: +1 888 826 9702  
Tel.: 610-448-4500  
Fax: 610-448-2254

#### Japan:

Siemens-Asahi  
Medical Technologies Ltd.  
Takanawa Park Tower 14 F  
20-14, Higashi-Gotanda 3-chome  
Shinagawa-ku  
Tokyo 141-8644  
(+81) 354238489

#### Asia:

Siemens Medical Solutions  
Asia Pacific Headquarters  
c/o Siemens Pte Ltd.  
Block 28 Ayer Rajah Crescent No. 06-06  
Singapore 139959  
(+65) 6341 0990

Or contact your local Siemens  
sales representative

Siemens AG · Medical Solutions  
Magnetic Resonance  
Henkestr. 127, D-91052 Erlangen  
Germany  
Telephone: ++49 9131 84-0  
www.SiemensMedical.com

**MAGNETOM Trio is the system of choice  
for true 3T whole-body applications.**

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases. The required features should therefore be specified in each individual case at the time of closing the contract.

Siemens reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens Sales representative for the most current information.

Note:  
Original images always lose a certain amount of detail when reproduced.

Siemens **Medical  
Solutions** that help

Order-No. A91100-M2220-B323-25-7600  
Printed in Germany  
CC 64323 SB 05035.