

UROSKOP Access

Top 10 Arguments

Customer Requirements

UROSKOP Access Benefits

■ Positioning of the patient only once and treatment of both kidneys without moving the patient	■ The first floor-mounted urodiagnostic X-ray table with truly symmetrical patient access from 4 sides of the table
■ Imaging of the whole urinary tract without moving the patient	■ Motorized movement of X-ray tube and Image Intensifier. Resting patient.
■ Optimized workflow (in terms of room planning, organizational flexibility, faster diagnostic routines and reduced consumption of anesthetics)	■ The patient is only positioned once and the anesthetist remains in one and the same clearly defined zone of the table, irrespective of what diagnostic procedure is being performed. Symmetrical patient access.
■ Digital real-time images of highest quality	■ Ideally matched imaging components with high resolution digital 1 K matrix yielding sharp, crystal-clear, high contrast images
■ Ergonomic image viewing	■ Two digital flat-screen monitors on flexible monitor arm for ideal viewing of images from any side of the table
■ Integration of endoscopy and urodynamics	■ Direct connection of endoscopy equipment and ultrasound/urodynamic systems via interface. UROSKOP monitors display endoscopy sequences in real time
■ Easy handling and full compatibility to the clinical network based on DICOM 3.0	■ Unparalleled integrated <i>syngo</i> look & feel operation, image post processing possibilities and MPPS. Full networking through DICOM 3.0 functionality.
■ Reliable table for accommodation of also very obese patients	■ Designed for patients up to 450 lbs (205 kg)
■ Ergonomic Trendelenburg tilt	■ Maintaining a constant working height in tilting around the perineal end of the table top
■ Maximum patient comfort	■ Table tilt of +/- 90° for images in the standing or sitting position. Seat height for micturition 19.7 inches (50 cm) only, like a normal toilet, very comfortable for the patient