

**ACUSON S2000™ Ultrasound System
Transducer Compatibility Matrix**

	ACUSON Sequoia System	ACUSON Antares, X300, X500, CV70 Systems & Sonoline G60 S System	ACUSON S2000 System	Biopsy Guide Information
18L6 HD			√	18L6 HD Needle Guide Set
14L5			√	SG-4 Needle Guide Set
14L5 SP			√	
9L4	√		√	SG-3 Needle Guide Set
4V1	√		√	SG-2 Needle Guide Set
4V1c	√		√	
6C2	√		√	SG-5 Needle Guide Set
4C1	√		√	CH4-1 Needle Guide Set
8V3	√		√	
4P1			√	
10V4	√		√	
EV8C4	√		√	Needle Guide available in future release
EC9-4			√	EC-1 Needle Guide Set, EC9-4 Disp. Needle Guide Set
7CF2			√	C7F2/C5F1 Needle Guide Set
9EVF4			√	EV9F4 Needle Guide Set
CW2		√	√	
CW5		√	√	
V5Ms	√	√	√	