



AX6ACOMNETCA - ACOM NET Supervisor

Course Number/Designator

AX6ACOMNETCA

Course Description

The ACOM NET Supervisor course is a two-day class covering the daily maintenance and operation of the ACOM NET System. Lectures include such topics as block diagram overview of the ACOM NET system, clinical operation, supervisor responsibilities check sheet, and supervisor console log on with details of all supervisor functions for the server, juke box ACOM M and ACOM P.C. clients. Lab sessions provide an opportunity for "hands-on" training in system operation, software configurations, database management, supervisor daily maintenance checks and media handling.

Upon completion of this course participants will be able to:

- Understand the operations of the ACOM NET functions and image data flow in a clinical environment.
- Perform daily checks through the supervisor computer.
- Change and maintain database backup tapes.
- Analysis of patient image database archive task completion.
- Recognize system deficiencies and utilize guidelines available through the supervisor web tool to take action for service reporting.
- Add or removal of MO disk media to juke box.
- Review work list tasks.
- Network deletion of Patient image data.

Audience

The ACOM NET Supervisor course is designed for the hospital personnel with the responsibilities of daily maintenance checks and database review.

Course/Skill Prerequisites

Experience in operation and daily use of the ACOM NET system and Cardiac Special Procedures lab workflow.

Course Information

- Classes begin at 8:30 AM each day and end at 5:00 PM (noon on last day).
- Tuition:** US \$2000 CE Credits: 10.5 hours, Category A
- Registration:** Fax a registration form (available online) to 1-919-468-7401 or call 1-800-888-SIEM (7436) opt 4.