

MammoTest Flx

Mobile Breast Biopsy System from Cintron Medical



Excellent imaging quality, maximum patient access, superior patient comfort, and now the **flexibility** to work where you want.

MammoTest Flx offers all the benefits of standard MammoTest®, plus:

- Mobile platforms for the table and generator allow the unit to be easily moved within a room or across locations
- Cable connections moved for more accessibility and quick disconnect added for convenience
- Table-control handles under the tabletop, decreasing table width and eliminating interference for the patient

Standard features include:

- Largest breast aperture in the industry at 11 inches/28 cm, ideal for arm-through procedures
- Rotating gantry and lateral arm to provide 360° patient and lesion access, especially suited for thin, compressed breasts and superficial lesions
- Flat detector, CCD camera (1024 x 1024 matrix), and CsI scintillator to allow for visualization of even the smallest microcalcifications

MammoTest Flx
Mobile Digital Prone Table Biopsy Solution

www.usa.siemens.com/healthcare

SIEMENS

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features which do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Local Contact Information

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway
Malvern, PA 19355-1406
USA
Phone: +1-888-826-9702
www.usa.siemens.com/healthcare

Global Business Unit

Siemens AG
Medical Solutions
Special Systems
Henkestr. 127
DE-91052 Erlangen
Germany
Phone: +49 9131 84-4746
www.siemens.com/healthcare

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
www.siemens.com/healthcare

Legal Manufacturer

Siemens AG
Wittelsbacherplatz 2
DE-80333 Muenchen
Germany

for the products

Cintron Medical Corporation
1130 W 124th Avenue, Suite 500
Westminster, CO 80234

Distributed by

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway
Malvern, PA 19355-1406
USA

www.siemens.com/healthcare